Expt. 25 Franck-Hertz experiment

I. Purpose

Verify the existence of atomic energy level through the Franck-Hertz experiment.

II. Principle

When an atom makes a transition from one stable state to another stable state under external actions, it absorbs or emits electromagnetic waves of a certain frequency, ν , given by:

$$h\nu = E_n - E_m$$

Where E_n and E_m are the n th and m th excited state energies of the atom respectively and h is Planck's constant.

The diagram of the experiment is shown as Figure 25 - 1. The Franck-Hertz Tube is a kind of tetrode, filled with argon gas. This experiment verifies the existence of atomic energy levels by investigating slow electrons colliding with atoms of argon. The hot cathode emits electrons and after being accelerated by the accelerating voltage U_{G_2K} between the cathode K and the second grid G_2 , they cross the grid G_2 . The decelerating voltage U_{G_2P} is applied between the plate and second grid G_2 . If the electron has enough energy, it can overcome U_{G_2P} to reach the plate, forming the plate current I_p . The main purpose of this experiment is to observe the variation of the plate current with accelerating voltage.


Figure 25 - 1 Schematic Diagram for Franck-Hertz Experiment

On gradually increasing the voltage U_{G_2K} , if the atomic energy level does exist, then the curve

$U_{G_2K} - I_P$ should be of form depicted in Figure 25 - 2. This curve shows the energy exchange between the argon atoms and the electrons in the space of $K - G_2$. The difference between two U_{G_2K} values corresponding to the two neighbouring minimum values of I_P equals the first excitation potential U_0 of the argon atom. When $U_{G_2K} = nU_0$, the plate current will decrease to a minimum value.


Figure 25 - 2 V - I Characteristic Curve in F - H Experiment

III. Experimental

1. Preheating

Before the start of experiment switch the measurement range to 10^{-6} , turning all potentiometers fully anticlockwise, to make the voltage to decrease to the minimum value. Then switch it on, setting the voltage selecting switch on the V_{G_2K} position and adjust (usually about 1.5V), and then switch to the V_{G_2P} position and adjust (usually about 7.5V). Begin observations after the F - H tube has preheated for three minutes.

2. Oscilloscope mode

① Switch the oscilloscope to the X - Y mode. Set the X axis amplifying ratio knob V/div to 0.2 V/div. The Y axis amplifying ratio V/div can be set to 20 mV/div or 50 mV/div. Then put the knobs to the CAL position (calibration), so that the actual amplifying ratios of the X axis and Y axis accurately equal the values indicated on the knobs. The "AC \perp DC" selecting switches for both X axis and Y axis should be in the DC position.

② Connect the output terminal of the F - H experiment instrument to corresponding input terminal of the oscilloscope. Adjust the gain control knob on the rear panel, making the horizontal scanning line in the screen cover 10 divisions, equivalent to 100 V.

③ Set the voltage selecting switch to the V_F position and adjust slowly. Once some wave crests are found with fast increasing I_A slightly decrease the value of V_F , till the wave form is stable.

④ Read 6 (or 5) voltage values corresponding to separate wave troughs. Then find the first excitation potential by the method of successive difference. With the Y axis of the oscilloscope standing for the plate current, you need to record the number of divisions to get your reading.

3. Manual mode

Observe the F - H curve, and adjust to find out the optimal filament voltage V_F . Switch the

voltage selecting switch to the position of "Accelerating Voltage", and adjust the accelerating voltage knob slowly. When the first peak value appears, record I_A and V_{C2K} , and then measure two nearby points of $V_{C2K} \pm 2V$; then, increase V_{C2K} , and when the first valley value appears, record I_A and V_{C2K} , then measure two nearby points of $V_{C2K} \pm 2V$. Totally measure 6 (or 5) peak values and 6 (or 5) valley values. Record the current value and the voltage value on the experiment instrument.

From the F - H experimental curve determine the first excitation potential.

4. Wave form of the accelerating voltage

Observe the wave form of the accelerating voltage on the second grid of the F - H tube, measuring its amplitude and frequency. Notice that the oscilloscope should be in the Y - t mode.

Notes: after the experiment is finished, turn the filament voltage V_F to the minimum (anti-clockwise), and then switch it off.

IV. Questions

1. In the condition of decelerating voltage $V_{C2P} = 0$, can you record a fluctuation of I_P ?

2. Consider the third wave trough of the F - H curve. At what location in the F - H tube does the relevant electron argon inelastic collision happen?

3. Based on the speed of the electrons emitted from the cathode, explain the form of the I_P wave crest? In the case of the initial velocity of all electrons being zero, would I_P decrease vertically? Draw the F - H curve for this situation?

4. For the manual mode and oscilloscope mode detail separately the characteristics of the voltage change on the second grid.

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