Expt. 1 Measurement of basic parameters of mechanics

I. Purpose

- 1) Master the use of conventional measuring instruments by measuring the density of regular- and irregular-shaped objects;
- 2) Measure the local acceleration due to gravity using the adjustable pendulum.

II. Apparatus

Vernier calipers, micrometer screw gauge, physical balance, adjustable pendulum, millisecond meter, thermometer, digital display caliper, pycnometer.

III. Principle

1. Measurement of the density of a solid substance (regular shape)

Assume that the mass of a uniform substance is m and the volume is V , so that the density is $\rho = m/V$. For a regular-shaped substance, m can be measured by using the physical balance, and V can be obtained indirectly by calculation of direct measurement with appropriate length measuring instruments, e. g., vernier calipers or micrometer screw gauge.

2. Measurement of the density of a substance (irregular shape)

For the irregular-shaped substance, V can be determined indirectly by Archimedes' principle. If the weight of the substance is $W_1 = m_1g$ in the air and $W_2 = m_2g$ under the liquid and we assume the density of the liquid is ρ' then the density of the measured substance should be

$$\rho = \frac{m_1}{m_1 - m_2} \rho' \quad (1)$$

When we use water as the medium and measure an irregular-shaped substance of density less than 1, we first determine the actual mass of the object and hang a weight underneath it, then immerse the weight completely to get the apparent mass m_2 of the system, as shown in Fig. 1. We then immerse the system completely to get the apparent mass m_3 . The density of the object is given by

$$\rho = \frac{m_1}{m_2 - m_3} \rho' \quad (2)$$

3. Measurement of the density of granular matter

Using pycnometer and distilled water the density of granular materials can be measured. The principle is still $\rho = m/V$, where the volume V needs


Figure 1 Measurement of the less-than-1 density material

to be obtained indirectly. There is a capillary in the cap of pycnometer. If the pycnometer is filled with water, then put on the cap, the excessive water will overflow. So the volume of matter in the pycnometer is always constant. In consequence, determine the mass of the granular matter measured m_1 , the mass of pycnometer filled with water ($\rho_0 = 1$) m_2 , and the mass the pycnometer filled with water and granular matter m_3 . Then the mass of water overflowed from the pycnometer is $m_1 + m_2 - m_3$, and the volume of the overflowed water is equal to that of m_1 . So the density of the measured granular matter is:

$$\rho = \frac{m_1}{m_1 + m_2 - m_3} \rho_0 \quad (3)$$

4. The measurement of the acceleration due to gravity g

Acceleration of gravity g is an important physical parameter in actual scientific research and manufacturing work. In this experiment, the measurement of the acceleration due to gravity can be done quickly and exactly by using the adjustable pendulum of which the exterior is shown in Fig. 2.

According to the formula of vibration cycle of pendulum $T = 2\pi\sqrt{\frac{L}{g}}$, when the pendulum lengths are L_i and L_j respectively, we can deduce:

$$g = \frac{4\pi^2(L_i - L_j)}{T_i^2 - T_j^2} \quad (1-4)$$

IV. Experimental

1. Measure the density of a solid using the physical balance

- (1) Adjust the horizontal bolt of the balance; make sure the air bubble is in the centre of the spirit level and the base plate of the balance is horizontal.
- (2) Zero the vernier to the left side of the beam. If the pointer is not central you should adjust the balance bolt and repeat this process, until the index pointer is at the centre of the dial.
- (3) Determine the density of the regular-shaped substance: put the object ($\rho > 1$), such as aluminum pieces or a steel ball, on the left scale pan and weights into the right pan. When the balance balances, read the number on the weights and vernier to give the mass of the object m_1 . Determine the volume of the aluminum pieces using vernier calipers, and the volume of the steel ball using the micrometer screw gauge, using the density definition formula calculate the density of two objects ρ also calculate the uncertainty range of the density, using the definition of the uncertainty for an indirect measurement.

- (4) Determine the density of an irregular-shaped substance ($\rho < 1$):

By the process described above determine the mass of the object m_1 ; place a beaker filled with water on the left scale pan, immerse the weight linked with the object completely into the water, read the apparent mass m_2 and immerse the system of weight and object completely into the water and read the apparent mass m_3 . From equation (3), determine the density, and work out the uncertainty of the result.

2. Measure the density of granular matter

based on the Archimedes principle, design procedures and measure the density of granular material using pycnometer and distilled water. Note: wipe off the water on the pycnometer each time.

3. Measure the local acceleration due to gravity by using the adjustable pendulum

In this experiment, the limit of intrinsic error of the ruler, which is used in measuring the pendulum length, is 0.02 mm, and that of the millisecond meter is 0.1 ms.

(1) The pendulum length differences $\Delta L = L_i - L_j$ should be chosen as 10 cm, 30 cm and 50 cm, in order to determine the effect of ΔL on the final precision of your measured value of g .

(2) Measure the period of the vibration of the simple pendulum, by the spreading method, which is: measuring 10 periods to give a good average for 1 period.

(3) Using equation (4) determine the local acceleration due to gravity " g " calculate the uncertainty range due to the precision of the apparatus.

Notes:

(1) Make sure that angle of swing of the simple pendulum is less than 5° and that it swings in a single plane.

(2) When changing the pendulum length, relax the screw of the line clamp. Tighten it again when making a measurement.

V. Questions

1. You are given a piece of alloy made of gold and copper. Describe how you could determine the ratio of gold to copper W_{Au}/W_{Cu} using the physical balance. Assume that the density of gold ρ_{Au} , the density of copper ρ_{Cu} and the density of water are known.

2. What kind of method could be used to increase the precision of the measurement of g , in this experiment?

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