## Expt. 23 Photoelectric effect and measurement of Planck constant

### I. Purpose

- (1) By using the photoelectric effect study the quantum nature of light.
- (2) To verify the Einstein equation and determine the Planck constant.
- (3) To obtain the voltage-current characteristics of a photoelectric tube.

### II. Apparatus

Mercury light and power supply, a range of light filters, diaphragm, photoelectric tube, bench, the Planck constant measuring instrument, connecting leads.

### III. Principle

When light with certain frequencies illuminates a metal electrons are emitted from the surface of that metal. This phenomenon is known as photoelectric effect. The kinetic energy of the emitted electrons does not depend on the intensity of light but does depend on the frequency of the incident light. This relationship can be written as:  $h\nu = E_{\text{kinetic}} + W$ , which is called the Einstein equation for the photoelectric effect.

In the Einstein equation, “ $h$ ”, the Planck constant is a very important physical quantity. In history, people have used many methods for its determination. Here, we will indirectly determine “ $h$ ” through the photoelectric effect.

As shown in Figure 23 - 1, A is the anode of the photoelectric tube, C is the cathode. G is a galvanometer, V is a digital voltmeter, R is a rheostat. By adjusting R, the voltage  $U$  between A and C can be changed continuously from negative to positive. If we measure the photoelectric currents at different voltages, then we can draw the voltage-current characteristics of the tube (see Figure 23 - 2).


Figure 23 - 1 Circuit Diagram of the Photoelectric Effect Apparatus


Figure 23 - 2 V - A Characteristics of the Photoelectric Tube

According to the conservation of energy, we have

$$h\nu = \frac{1}{2}mv_{\max}^2 + h\nu_c \quad (23-1)$$

This formula shows the basic physics of the photoelectric effect. In the Einstein equation,  $h\nu$  is the energy of incident photon,  $m$  is the mass of the electron,  $v_{\max}$  is the maximum velocity of the emitted electron,  $\nu_c$  is the threshold frequency of the photoelectric effect, and  $h\nu_c$  is the work function of the metal. The threshold frequency varies from metal to metal. If we apply an increasing reverse voltage on electrode A (anode) until the emitted electrons which have maximum kinetic energy are prevented from reaching the anode then:  $eU_c = E_k = \frac{1}{2}mv_{\max}^2$ , where  $U_c$  is called the stopping voltage.

Therefore  $U_c = h(\nu - \nu_c)/e \quad (23-2)$

If we measure values of  $U_c$  for various situations experimentally then  $h$  can be determined.

#### IV. Experimental

The photoelectric effect apparatus is shown in Figure 23-3.


Figure 23-3 Photoelectric Effect Apparatus

##### 1. Preparing

Turn on the supply, and allow it to warm up for 20 minutes. Cover the light sources with the hood. Adjust the distance between the photoelectric tube and the mercury light to about 40cm. Connect the voltage input terminals by connecting leads of the same colour.

##### 2. Measurement of Planck constant (tube A)

Considering the stability of the experiment, we choose the zero-current method. We record the absolute value of the voltage  $U$  as the stopping voltage  $U_c$  at a certain frequency of illumination when the measured current is zero.

Set the voltage button on the  $-2 \sim +2V$  level, and the current button on the  $10^{-12}A$  level. Remove the current input leads and reconnect after zeroing in the instrument.

Use the diaphragm with a 4mm diameter and a light filter on the entrance hole of photoelectric tube. The tube and the mercury light are  $L=400mm$  apart. Turn the voltage-adjusting knob anticlockwise; observe the variation of current while decreasing the voltage between A and C from a positive value to  $-2V$ , and record the values of stopping voltage (as  $U_c$ ) when the photocurrent reaches zero.

Measure and construct  $U - V$  curve. Compare the calculated  $h$  value with the accepted value of the Planck constant  $h_0$ .

Work out the relative error  $E = (h - h_0) / h_0$ , where  $h_0 = 6.626 \times 10^{-34} \text{ J} \cdot \text{s}$ . So far, we have measured the Planck constant  $h$ .

### 3. Measuring the $U - I$ characteristic curve of the photoelectric tube (tube B)

Turn on the supply; set the voltage on the  $-2 \sim +30\text{V}$  level and the current range on the  $10^{-9}\text{A}$  level. Use the diaphragm with a 4mm diameter and a light filter on the entrance hole of photoelectric tube. The tube and mercury light are  $L=400\text{mm}$  apart. Do the same zeroing as mentioned in procedure 2. Measure the required data and fill in the following table; construct the  $U - I$  curve.

### 4. Verifying the relationship of the saturation current and the intensity of the incident light

Set  $U=30\text{V}$ . Do the same as described in procedure 2. Measure the required data and construct the  $I - P$  curves.

**Notes:** You must cover the output hole of the light source with the hood while changing the light filter. After you finish the entire experiment, please cover the input hole of the photoelectric tube. This avoids high intensity light illuminating the cathode and shortening the lifetime of the tube.

## V. Questions

1. How did you determine the stopping voltage? Explain your method.

2. Can the light filter be inserted into the entrance of the diaphragm while doing the experiment?

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