Expt. 8 Hall effect and measurement of parameter

I. Purpose

- (1) Use the Hall Effect to investigate the performance of semiconductor materials.
- (2) Use the "Symmetry Measuring Method" to eliminate the influence of secondary effects.

II. Apparatus

Hall Effect experiment instrument, Hall Effect experiment combiner, Hall Effect magnetism measuring instrument, teslameter, multimeter, etc.

III. Principle

1. Hall Effect

If an electric current flows through a conductor in a magnetic field, the magnetic field exerts a transverse force on moving charges which tends to push them to one side of the conductor. A build up of charges at sides of the conductor results in additional transverse electric field in a direction perpendicular to the direction of magnetic field. This phenomenon is called Hall Effect.

As shown in Figure 8-1, we apply a current I_s in the y direction of the semiconductor sample and we apply a magnetic field B in z direction. This causes electric charge to accumulate on the x axis, generating a Hall voltage U_H . The electric field E_H caused by the electric potential difference is called the Hall field. If we assume that the carrier concentration is n , then the relation between the Hall voltage U_H , the magnetic induction density B and other quantities such as carrier concentration, n , is


Figure 8-1 Principle of the Hall Effect

$$U_H = \frac{I_S B}{ned} = R_H \frac{I_S B}{d} = K_H I_S B \quad (8-1)$$

Where: $R_H = 1/ne$, R_H is called the Hall coefficient. This is an important parameter which reflects magnitude of the Hall Effect in the material. $K_H = R_H/d = 1/ned$, where K_H is Hall element sensitivity, unit in $V/(A \cdot T)$.

The positive or negative sign of the Hall voltage (or value for R_H) and the direction of magnetic field determine the conduction type of the sample. Semiconductor materials can be divided into P type (hole type) and N type (electron type). If the measured $U_H > 0$, then R_H is positive and the sample is P type. If $U_H < 0$, then R_H is negative and the sample is N type.

The Hall element sensitivity K_H is a constant, in general the bigger the better. If we measure the values of the Hall current I_S and corresponding Hall voltages U_H we can determine the magnitude of the magnetic induction density B ; this is the principle of the measurement of magnetic field using the Hall Effect.

The conductivity σ of a semiconductor is given by:

$$\sigma = \frac{I_S L}{U_{CA} S} \quad (8-2)$$

Carrier mobility μ : the relationship between conductivity σ , carrier concentration n and mobility μ is as follows:

$$\mu = |R_H| \sigma \quad (8-3)$$

2. Elimination of secondary effects in the Hall Effect

The Hall Effect is accompanied by multiple secondary effects which create voltages superimposed on the Hall voltage, causing errors to the measurement of the Hall Effect. These secondary effects can all be decreased or eliminated by changing the direction of the operating current or magnetic field. We change the current direction and the magnetic field direction in turn and take measurements under four conditions: $(+B, +I_S)$, $(-B, +I_S)$, $(+B, -I_S)$ and $(-B, -I_S)$. We take the average value of each "Hall voltage" to give:

$$U_H = \frac{1}{4} (|U_1| + |U_2| + |U_3| + |U_4|) \quad (8-4)$$

3. Coupling degree of a pair of coaxial coils

When the span " a " of a pair of coaxial coils equals the radius R of each coil, they form the so called "Helmholtz coils". These coils have a uniform distribution of magnetic field intensity H on the axis but when $a \neq R$, the H on the axis will be non-uniform. The coupling degree of the two coils (for example under-coupling, over-coupling) can be measured by the Hall Effect.

IV. Experimental

(1) Measure the electromagnet excitation curve, i. e. $I_M - B$ Curve, by Teslameter.

The Teslameter is the instrument which uses the principle of the Hall Effect to measure magnetic induction density, B . Measure the magnetic induction density B at the centre of the electromagnet

for the magnet current values I_M (0.100 – 0.800 every 0.100). Draw the $I_M - B$ curve.

Notes: The Hall probe is made of extremely thin semiconductor material, which is very fragile. Be very careful when using it! Protect it by the casing tube whenever it is not in use.

(2) Fixing I_M , measure $U_H - I_S$ curve.

Take $I_M = 0.50A$. Adjust I_S , the sample current, making I_S at certain value, such as $I_S = 1.00mA$, then using the method to eliminate the secondary effects; get U_H values in the range of $I_S = 1.00 - 10.00mA$. Measure corresponding U_H values for each increase in 1.00mA, and draw the $U_H - I_S$ relation curve.

(3) Fixing I_S , measure $U_H - I_M$ curve.

Make $I_S = 5mA$. Adjust I_M , making $I_M = 0.100A, 0.200A, \dots, 0.800A$. Use the method to eliminate secondary effects to get corresponding U_H value and draw the $U_H - I_M$ relation curve.

(4) Under zero magnetic field ($B = 0$), $I_S = 0.20mA$, measure the value of U_{CA} .

(5) Based on the formula, calculate the Hall sensitivity K_H , the Hall coefficient R_H and the carrier concentration n , finally determine $\overline{K_H}$, $\overline{R_H}$, \overline{n} .

(6) Use the formula to calculate the values of conductivity and carrier mobility. For the Hall plate $b = 4mm$, $L = 3mm$, $d = 0.5mm$.

(7) Determine the conduction type of the sample from the sign of the Hall voltage.

Whether it is N type or P type?

(8) Observe the coupling degree of a pair of coaxial coils (see Figure 8 - 2).


Figure 8 - 2 Relation of Distribution of Magnetic Field and Coil Spacing on the Axis

First of all connect the two coaxial coils in series, then adjust the exciting current of the Hall Effect combiner to 0.80A. Set the current of Hall Effect magnetism measuring instrument to 8mA. Alter the distance "a" between the two coaxial coils; by making, in turn, $a > R$, $a = R$, $a < R$. Measure the values of "a" and R.

Move the probe bar of the Hall element along the axial direction between the two coils. Record the probe bar travelling distance X and the corresponding Hall voltage value. Take more than 10 different position points for each of the three coupling states separately, and design the data table by yourself. Draw the $U_H = f(x)$ relation curve of the coaxial coils in the three coupling states.

Notes: The Hall plate only uses a very small current. NEVER mistakenly connect it to the exciting current of the Hall Effect Combiner!

V. Questions

1. What precautions should you take when using a teslameter to measure magnetic field?
2. What requirements are there on a material used to manufacture a Hall plate?

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