Expt. 10 Interference of light

I. Purpose

(1) Observe two kinds of interference phenomena: wedge interference and Newton Rings interference.

(2) Use the principle of wedge interference to measure the diameter of glass fiber; measure the radius of curvature of a lens by Newton Rings.

II. Apparatus

Measuring microscope, sodium vapour lamp, Newton Rings assembly, optically flat glass.

III. Principle

When two trains of monochromatic light with constant oscillation direction, frequency and phase difference are superimposed, the resulting oscillation is reinforced in some places and weakened in others. This is called the interference of light.

Interference is widely used in scientific research and engineering. For instance, interference of light can be used to accurately measure length and its variation, detect the quality of optical element surfaces, determine spectral wavelength and their fine structure, etc.

1. Wedge interference

Place the glass fiber between two pieces of flat glass, parallel to their intersecting edges as shown in Figure 10 - 1. When a monochromatic parallel beam of light arrives perpendicular ($i=0$) at the air wedge ($n=1$) formed by the two pieces of flat glass, the two beams of reflected light a' and b' at the wedge tip C cause interference, forming fringes of interlaced light and shade. According to the thin film interference formulae, here


Figure 10 - 1 Principle of Wedge Interference

$$\begin{cases} \delta = 2e + \frac{\lambda}{2} = 2k \frac{\lambda}{2} & k = 1, 2, 3, \dots \text{Bright fringe} \\ \delta = 2e + \frac{\lambda}{2} = (2k + 1) \frac{\lambda}{2} & k = 0, 1, 2, \dots \text{Dark fringe} \end{cases}$$

The thicknesses of air at various points of the same interference fringe are equal, so these are called interference fringes of equal thickness. Where the two pieces of glass touch, $e = 0$, the optical path difference of the two beams of light $d = \lambda/2$, so there should be dark fringe. Take the air thickness of the k_1^{th} dark fringe as e_1 , the air thickness of the $(k_1 + \Delta k)^{\text{th}}$ dark fringe as e_2 , the horizontal distance of the two fringes as x , and the air thickness difference between the dark fringe as Δd , then we have:

$$\Delta d = e_2 - e_1 = \Delta k \frac{\lambda}{2}$$

$$\tan \alpha = \Delta d / x = \Delta k \lambda / 2x$$

Glass fiber diameter

$$D = 1 \cdot \tan \alpha = 1 \lambda \Delta k / 2x.$$

2. Newton Rings

Newton Rings are shown in Figure 10-2.

When the parallel light beam arrives at the air wedge formed by the lower surface of a lens of very large curvature radius and the upper surface of the optically flat glass we get interference of light. The interference fringes are concentric circles around the point of contact of the lens and the glass due to equal thickness interference. They are Newton Rings.

Using the same interference principle as wedge interference, through analysis we can find that the relationship between the radius r_k of the k^{th} dark fringe of the Newton Rings and the radius R of the lens curved surface is:

$$R = \frac{r_k^2}{k\lambda} \quad (k = 0, 1, 2, \dots)$$

Notes: the centre of Newton Rings is not theoretically a dark point, but a dark spot, which means we cannot determine the exact centre and order k of the Newton Rings. So we can take k as any values of m and n (for example, $m = 11$, $n = 1$) so that

$$\begin{cases} r_m^2 = mR\lambda \\ r_n^2 = nR\lambda \end{cases}$$

This gives:

$$R = \frac{r_m^2 - r_n^2}{(m - n) \lambda}$$

The wavelength of the sodium vapor lamp used in the experiment is $\lambda = 589.3 \text{ nm}$.

IV. Experimental

Under the microscope of the experiment system, there is a semi-reflecting mirror that can reflect the parallel light to the microscope workbench. By turning drum wheels workbench can be moved along X axis and Y axis respectively. The drum wheel can be read to 0.01mm.


Figure 10-2 Setup of Newton Rings

1. Observation of wedge interference

(1) Place the glass plate on the microscope workbench, and clip a piece of glass fiber between the two pieces of glass.

Notes: ① make the glass fiber parallel to the edge; ② away from the wedge tip; ③ the direction of the wedge tip should be perpendicular to the moving direction of the workbench.

(2) Focus the microscope, until you can see clear interference fringes.

(3) Make the intersection point of the cross hair nearly at the end of the wedge. Record the position of a certain dark fringe, then count 30 dark fringes (take Δk as 30), and measure the distance X of the 30 dark fringes. Repeat this experiment for five times continuously counting in the same direction, totally counting 5×30 fringes. This avoids the backlash effect of back and forth motion of the microscope.

(4) Measure the distance L from the wedge tip to the glass fiber.

2. Observation of Newton Rings

(1) Place the Newton Rings assembly on the workbench, turn on the sodium vapour lamp and set the semi-reflecting mirror, so as to see the visual field through the eyepiece as bright as possible.

(2) Adjust the eyepiece, so you can see the cross hair clearly, and then focus on the Newton Rings.

(3) Adjust the drum wheel, so that the cross hair goes through the central spot of the interference pattern.

(4) Turn the drum wheel, to align the intersection point of the cross hair on the first ring surrounding the Newton Rings dark spot. Then turn the drum wheel again, and count to the 11th ring moving toward right. To eliminate the effect of backlash (lost motion), you should move the cross hair further to the right, and then return to the position of the 11th ring. Record the position (X_{11}) here, and then move leftward to the 1st ring, record the position (X_1). Continue to move leftward to the 1st ring at the other side of the central spot (recorded as X_1'), and then continue to move leftward to the 11th ring (recorded as X_{11}'). In this way you have measured 4 positions and from $(X_{11} - X_{11}')/2$ and $(X_1 - X_1')/2$ you can get r_{11} and r_1 . Repeat above steps six times.

V. Questions

1. In the Newton Rings experiment, what is the difference between the interference fringes formed by reflected light and transmitted light?

2. Interference fringes generated in the wedge interference experiment are not totally parallel to the edge. Explain what causes this?

3. When measuring the plano-convex lens radius of curvature R using Newton Rings does it matter if you do not cross the centre of the circular interference fringes while measuring r_m and r_n ?

Can you still use the formula $R = \frac{r_m^2 - r_n^2}{(m - n)\lambda}$?

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